testo 206  
pH-/Temperatur-Messgerät

Bedienungsanleitung	de
Instruction manual	en
Manuale di istruzioni	it
Руководство по эксплуатации	ru

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## General Information

Please read this document through carefully and familiarise yourself with the operation of the product before putting it to use. Keep this documentation to hand so that you can refer to it when necessary.

### Symbols and what they mean

Symbol	Meaning	Remarks
 <b>Warning!</b>	Warning text: Warning! Serious physical injuries could occur carefully and take the if the precautionary measures specified are not taken.	Please read warning precautionary measures specified.
 <b>Caution!</b>	Warning text: Caution! Light physical injuries could occur if the precautionary measures if the precautionary measures specified are not taken.	Please read warning carefully and take the precautionary measures specified.
Note 		Pay particular attention to Notes.
	Button name	Press button.
Text, 	Display content	Text or symbol is shown on the display.

# Contents

General Information .....	20
Contents .....	21
1. Safety Information.....	22
2. Intended Purpose .....	23
3. Product Description.....	24
3.1 Display and control elements.....	24
3.2 Probe/BNC modules .....	24
3.3 Power.....	25
3.4 TopSafe .....	25
3.5 Storage cap .....	25
3.6 Wall/Transport bracket .....	25
3.7 Replacing the storage cap.....	25
4. Initial Operation.....	26
4.1 Insert battery .....	26
4.2 Connect external probe (BNC module pH3 only).....	26
5. Operation .....	26
5.1 Switching on/off .....	26
5.2 Setting instrument .....	26
5.3 Opening setting mode .....	27
5.4 Measuring .....	29
5.5 Calibrating instrument.....	31
6. Service and Maintenance .....	32
6.1 Checking electrolyte gel.....	32
6.2 Cleaning housing.....	32
6.3 Cleaning probe.....	32
6.4 Disinfecting the probe .....	33
6.5 Changing the module .....	33
6.6 Changing battery.....	33
7. Questions and Answers.....	34
8. Technical data .....	35
9. Accessories and Spare Parts.....	36

## 1. Safety Information

### **Avoid electrical hazards:**

- ▶ Do not measure on or near live parts!

### **Preserving product safety/warranty claims:**

- ▶ Operate the instrument properly and according to its intended purpose and within the parameters specified. Do not use force.
- ▶ Do not store together with solvents (e.g. acetone).
- ▶ Temperature data on probes/sensors refer only to the measurement area of the sensors. Do not expose handles and cables to temperatures greater than 70°C if they are not specifically designed for higher temperatures.
- ▶ Open the instrument only when this is expressly described in the documentation for maintenance purposes.
- ▶ Only the maintenance and service work described in the documentation should be carried out. Please adhere to the steps specified. For safety reasons, only original spare parts from Testo should be used.

### **Ensure correct disposal:**

- ▶ Disposal of defective rechargeable batteries and spent batteries at the collection points provided.
- ▶ Send the instrument directly to us at the end of its life cycle. We will ensure that it is disposed of in an environmentally friendly manner.

### **Protect the environment:**

The button cell used in the instrument contains 1,2-Dimethoxyethane (CAS 110-71-4). See EC Regulation No. 1907/2006 (REACH) Art. 33.

## 2. Intended Purpose

testo 206 is a practical instrument for spot check measurements of pH value and temperature.

It has different applications depending on the probe/BNC module used.

### testo 206 with pH1 immersion probe

Measurement of liquid substances in the following sectors:

- Food companies (e.g. fruit juices)
- Industry (e.g. refrigerant, electroplating, chip production, paints and varnish, print products)
- Chemistry (e.g. cleaning agents)
- Environmental protection (e.g. drinking water / waste water)
- Swimming pools, aquaria
- Agriculture
- Fish farming
- Pharmaceuticals and biotechnology

### testo 206 with pH2 penetration probe

Measurement of semi-solid substances in food production and processing: e.g. marmelades, marzipan, pastes, ready-to-serve salad, gelling agents, fruit, milk products, bakery and confectionery products. Laboratory measurements in companies involved in the processing of food.

### testo 206 with pH3 BNC module

The BNC socket is used to connect external pH probes. Ranges of application depend on the probe attached.

testo 206 is not suitable for diagnostic measurements in the medical sector.


The following components of the product are designed for continuous contact with foodstuffs in accordance with the regulation (EC) 1935/2004:

The measurement probe up to 1 cm before the probe handle or the plastic housing. If provided, the information about penetration depths in the instruction manual or the mark(s) on the measurement probes should be noted.

## 3. Product Description

### 3.1 Display and control elements

### 3.2 Probe/BNC modules

Immersion probe (pH1)

Penetration probe (pH2)

BNC module (pH3)

### 3.3 Power

Power is via a button cell (Typ CR2032, 3V; included with delivery).

### 3.4 TopSafe

TopSafe protects the instrument from moisture and mechanical strain (impact). We recommend you always use TopSafe.

! The IP68 protection class is achieved if the instrument is placed inside the Top-Safe and is closed.

### 3.5 Storage cap

The storage cap filled with electrolyte gel is used to store the probe between measurements.

The probe is only immediately ready to operate if stored in electrolyte gel.

If the probe has been out of the electrolyte gel for a longer period of time, it must be stored in the electrolyte gel for approx. 12 hours to regenerate.

The storage cap can also be attached to the wall/transport bracket.

### 3.6 Wall/Transport bracket

The wall/transport bracket with belt holder and attachment unit for the storage cap is used for the safe storage of the measurement instrument at a fixed point or during transport.

### 3.7 Replacing the storage cap

Remove the instrument from the wall/transport bracket. The storage cap is simply plugged on and can be pushed off the bracket by applying light pressure.

## 4. Initial Operation

### 4.1 Insert battery

1. Open the TopSafe sealing cap and remove instrument.
2. Open battery compartment on rear of instrument.
3. Insert button cell (Type CR2032, 3V). The (+) sign must be visible.
4. Close battery compartment.
5. Remove protection strip on the storage cap.

### 4.2 Connect external probe (BNC module pH3 only)

- ▶ Attach BNC plug from the external probe to the BNC socket and seal with the aid of a bayonet coupling.

## 5. Operation

### 5.1 Switching on/off

- ▶ Switch on instrument: .
  - All segments light up briefly and the instrument changes to the measurement mode.
- ▶ Switch off instrument: Keep  pressed.

### 5.2 Setting instrument

Overview of functions:

Function	Description	Einstellmöglichkeiten
Temperature unit	Sets unit	°C or °F
Auto Hold ( <b>AUTO HOLD</b> )	Automatically holds reading, as soon as it is stable*	<b>On</b> (switched on) or <b>OFF</b> (switched off)
Gradient/Offset	Displays gradient and offset values stored in the instrument (amount)	None (Information only)
Calibration method ( <b>CAL</b> )	Sets 1, 2 or 3 point calibration	<b>1P</b> , <b>2P</b> or <b>3P</b>
Calibration points ( <b>CAL pH</b> )	Sets calibration points	1P: <b>4</b> , <b>7</b> or <b>10</b> 2P: <b>4 7</b> or <b>7 10</b>
Auto Off ( <b>AUTO OFF</b> )	Instrument switches off automatically after 10 minutes if no button has been pressed	<b>On</b> (switched on) or <b>OFF</b> (switched off)
Beeper ( <b>bP</b> )	Warning sound (button pressed, stable reading reached with	<b>On</b> (switched on) or <b>OFF</b> (switched off)

\* Modification less than 0.02pH in 20s

The following functions can be set:

## 5.3 Opening setting mode

To access setting mode, the instrument must be switched off.

1. Open setting mode: Hold down MODE + tap ON/HOLD.

### 5.3.1 Setting the temperature unit (°C or °F)

1. Open setting mode.
2. °C or °F appears.
3. Select temperature unit via CAL.
4. Press MODE to confirm selection
5. Via MODE, go through the settings until all segments light up briefly.
6. The instrument then switches to the measuring mode.

### 5.3.2 Switching Auto-Hold on or off

1. Open setting mode.
2. Via MODE, go through the settings until AUTO HOLD appears.
3. Select On or OFF via CAL.
4. Press MODE to confirm selection.
5. Via MODE, go through the settings until all segments light up briefly.
6. The instrument then switches to the measuring mode.

### 5.3.3 Displaying the gradient and offset values stored in the instrument

1. Open setting mode.
2. Via MODE, go through the settings until mV/pH (gradient and offset values) appear.
3. Via MODE, go through the settings until all segments light up briefly.
4. The instrument then switches to the measuring mode.

### 5.3.4 Selecting the calibration method (1P, 2P or 3P method)

1. Open setting mode.
2. Via MODE, go through the settings until CAL 1P, 2P OR 3P appears.
3. Select calibration method 1P, 2P OR 3P via CAL
4. Press MODE to confirm
  - 1P: can be selected via CAL (4), (7) or (10)
  - 2P: can be selected via CAL (4,7) or (7,10)
5. Press MODE to confirm selection
6. Via MODE, go through the settings until all segments light up briefly.
7. The instrument then switches to the measuring mode.

### 5.3.5 Switching Auto Off on and off

1. Open setting mode.
2. Via MODE, go through the settings until AUTO OFF appears.
3. Select On or OFF via CAL
4. Press MODE to confirm selection
5. Via MODE, go through the settings until all segments light up briefly.
6. The instrument then switches to the measuring mode.

### 5.3.6 Switching the beeper on/off

1. Open setting mode.
2. Via MODE, go through the settings until bP appears.
3. Select ON or OFF via CAL
4. Press MODE to confirm selection
5. Via MODE, go through the settings until all segments light up briefly.
6. The instrument then switches to the measuring mode.

## 5.4 Measuring

### Setting up the instrument

#### Warning!

Measurement tip made of glass, risk of breakage!

Risk of injury on account of glass parts which remain in the measurement medium.

- ▶ Check measurement tip of pH probe after each measurement for damage.

**!** If large quantities of the electrolyte gel are stuck to the probe when removed from the storage cap, it is a sign that the gel is spent.

- ▶ A new storage cap is needed.
- ▶ Clean off the pH probe before and after each measurement with low concentration soap water followed by tap water (water temperature should be below 40 °C).
- ▶ When using the BNC module, please note the application information included with the external probe.

After horizontal storage:

- ▶ Briefly shake the probe in order to release any gas bubbles which may have formed in the probe tip.

1. Carefully remove storage cap.
2. Switch on instrument: .

### Carry out measurement

- ▶ Immerse/penetrate probe in the medium to be measured.
- The measured pH and temperature values are displayed. The readings are updated twice a second.
  - ▶ Hold readings manually: .
  - ▶ Restart measurement: .
  - If Auto-Hold is switched on, AUTO HOLD flashes until the instrument has found a stable pH reading. The readings are then frozen (AUTO HOLD lights up). If no stable pH readings have been found within 300s, the measurement is stopped (⊖ and AUTO HOLD light up).
  - ▶ Restart measurement: .

**! Manual temperature compensation**

The manual temperature compensation only influences the reading if a pH probe without a temperature sensor is connected. The instrument can therefore be adjusted to the temperature of the measurement medium.

This function is only available with a connected BNC module (pH3) if a pH probe without temperature sensor is connected. The temperature can then be adapted to the temperature of the medium being measured.

1. Open manual temperature compensation mode: **MODE**.
  - ▶ Increase value: **CAL**. Keep button pressed to get through values quickly.
2. Change setting direction: **MODE**.
  - ▶ Lower value: **CAL**. Keep button pressed to get through values quickly.
3. Finish setting: **MODE**.
  - All segments light up briefly and the instrument will change to the measurement mode.

**Finish measurement**

1. Switch off instrument: Keep **ON/HOLD** pressed.
2. Clean off the pH probe with low concentration soap water followed by tap water (water temperature should be below 40 °C).
3. Insert probe in the storage cap.

**!** The probe tip must be immersed in the electrolyte gel.  
Keep electrolyte gel clean.

## 5.5 Calibrating instrument

! Please also adhere to the instructions supplied with the buffer solution (Testo buffer: see label).

! During calibration it is important, that the glass probe does not touch the synthetic material of the bottle. Preferably do not leave the instrument in the bottle, as variations in calibrations of up to  $\pm 0.4$  pH can arise.

! The instrument is switched on and is in the measurement mode.

1. Open calibration mode: .

  - testo 206-pH3 with pH probe without temperature sensor: The set temperature value for manual temperature compensation is displayed for 2s. The value must correspond to the temperature of the buffer solution.
  - Calibration point (4, 7 or 10) is shown and CAL flashes.

2. Skip calibration point: .

  - or-
  - Immerse probe in the buffer solution and start calibration: .

    - The instrument waits for a stable reading : AUTO flashes.
    - If a stable reading is available (change less than 0.02pH in 20s), the calibration point calibrates and the instrument changes to the next calibration point (if available) or to the gradient and offset value display.
      - ▶ Carry out calibration manually: .

  - 3. Repeat Step 2 for additional calibration points.
    - Once calibration is complete, the amount of the gradient and offset value is shown. If the amount of the gradient value is less than 50mV/pH or the amount of the offset value is greater than 60mV, the pH electrode is spent and must be replaced.
  - 4. Return to measurement view: .

## 6. Service and Maintenance

### 6.1 Checking electrolyte gel

- ▶ Check the electrolyte gel in the storage cap regularly for contamination and ensure it is filled to the correct level. Replace storage cap if necessary.

### 6.2 Cleaning housing

- ▶ Clean housing with a damp cloth (soap water) if dirty. Do not use abrasive cleaning agents or solutions! TopSafe can also be cleaned in the dishwasher.

### 6.3 Cleaning probe

---

Destruction of probe due to incorrect cleaning!

Danger of injury due to glass parts remaining in the measurement medium.

- ▶ Use only the cleaning substances stated.
- 

Depending on the type of contamination, the following cleaning substances are suitable:

- Fat: domestic dishwashing liquid
- Protein: pepsin

The use of warm water improves the cleaning effect.

1. Apply dishwashing liquid or pepsin to a cloth and gently wipe (do not rub, as this leads to static charging).
2. Rinse the probe with clear warm water.
3. To stabilize the probe, place in the storage solution for at least 1 hour (better 12 hours).
4. Re-calibrate the probe (see 5.5 Calibrating instrument).

## 6.4 Disinfecting the probe

The measuring instrument should not be exposed to disinfectant for an excessive period of time. Using the TopSafe protective case can reduce the risk of incompatibility.

1. Wet a paper towel with disinfectant.
2. Gently dab the probe with the disinfectant, avoid rubbing.
3. Ideally, regenerate the probe in electrolyte gel for 12 hours before carrying out another measurement. (Using alcoholic disinfectants removes water from the glass membrane)
4. After each disinfection, carry out calibration with the buffer solution. (see 5.5 Calibrating the instrument)

## 6.5 Changing the module

**!** The instrument must be newly calibrated if a module is changed (See 5.5 Calibrating instrument)!

Instrument must be switched off. Plug contacts in the instrument should not be touched!

1. Unscrew the screws at the back of the instrument.
2. Remove module and attach new module.

**!** There must be rubber sealing rings available for the screws.

- ▶ Check that the sealing rings are in their correct position.

3. Tighten screws.

## 6.6 Changing battery

1. Open battery compartment at the back of the instrument.
2. Remove spent button cell and insert new button cell (Type CR2032, 3V). The (+) symbol should be visible.
3. Close the battery compartment.

## 7. Questions and Answers

Question	Possible causes	Possible solution
Readings instable.	Static charge.  Air cushion from measurement electrode gets into measurement tip. pH electrode has dried out.	<ul style="list-style-type: none"> <li>▶ Rinse off pH electrode with tap water or low concentration soap water.</li> <li>▶ Shake pH electrode in a downward direction like with a fever thermometer.</li> <li>▶ Place pH electrode for several hours in water or diluted hydrochloric acid.</li> </ul>
 lights up.	Remaining capacity of battery < 10h.	<ul style="list-style-type: none"> <li>▶ Change battery (See 6.6 Changing battery 33)</li> </ul>
Instrument switches itself off.	Auto Off function is switched on.	<ul style="list-style-type: none"> <li>▶ Auto Off. (See 5.2 Setting instrument 26)</li> </ul>
<b>Er1</b> lights up.	Invalid gradient value of pH electrode.  pH electrode defective.	<ul style="list-style-type: none"> <li>▶ Recalibrate instrument, a new buffer solution could be used.</li> <li>▶ When using a probe without temperature sensor: check set temperature value.</li> <li>▶ Change probe.</li> </ul>
<b>Er2</b> lights up.	Invalid offset value of the pH electrode.  pH electrode defective.	<ul style="list-style-type: none"> <li>▶ Recalibrate instrument, a new buffer solution could be used.</li> <li>▶ Change probe.</li> </ul>
<b>Er3</b> lights up.	Invalid gradient value of pH electrode following 3 point calibration. pH electrode defective.	<ul style="list-style-type: none"> <li>▶ Recalibrate instrument, a new buffer solution could be used.</li> <li>▶ Change probe.</li> </ul>

If we have not answered your question, please contact your nearest distributor or Testo Customer Service point.

For contact data, see back of this document or web page [www.testo.com/service-contact](http://www.testo.com/service-contact).

## 8. Technical data

Typ	testo 206-pH1	testo 206-pH2	testo 206-pH3
Parameters		pH/°C	
Sensor		pH electrode/NTC	
Measurement range		0 ... 14 pH/±0 ... +60 °C (short-term to +80°C, max. 5min)	
Resolution		0.01 pH/0.1 °C	
Accuracy		±0.02 pH/±0.4 °C	
Temperature compensation	Automatic	Automatic	Depending on external probe
Probe	Probe module w/immersion probe	Probe module w/penetration probe	BNC module w/ connection socket
 Measuring rate		2/s	
Operating temperature		±0 ... +60 °C	
Storage temperatur		-20 ... +70 °C	
Power		1 x button cell, Type CR2032, 3V	
Battery life		Approx. 80h	
Housing		Instrument: ABS, TopSafe: PU	
Protection Class		With TopSafe: IP 68	
CE guideline		2014/30/UE	
Dimensions (L x B x H)		110 x 33 x 20 (without probe and Topsafe)	

## 9. Accessories and Spare Parts

Name	Item no.
Replacement pH probe for testo 206 pH 1 including gel storage cap	0650 2061
Replacement pH probe for testo 206 pH 2 including gel storage cap	0650 2062
Storage cap for testo 206 pH 1/pH 2 with KCl gel filling	0554 2067
pH universal plastic electrode without temperature sensor for testo 206 pH 3, including storage/wetting cap	0650 2063
pH universal plastic electrode with temperature sensor for testo 206 pH 3, including storage/wetting cap	0650 2064
pH glass electrode with temperature sensor for testo 206 pH 3, including storage/wetting cap	0650 1623
pH food electrode without temperature sensor for testo 206 pH 3, including storage/wetting cap	0650 0245
Electrolyte solution (50 ml) for storing pH electrodes in the storage/wetting cap	0554 2318
Replacement storage/wetting cap (50 ml) for pH electrodes	0554 0048
Replacement lithium button cell type CR 2032	0515 0028
pH buffer solution 4.01 in dosing bottle (250 ml), including DAkKS calibration certificate	0554 2061
pH buffer solution 7.00 in dosing bottle (250 ml), including DAkKS calibration certificate	0554 2063
ISO calibration certificate analysis for pH buffer solutions; calibration points 4 pH, 7 pH, 10 pH	0520 0007
ISO calibration certificate analysis at 3 pH values over the measuring range	0520 0037

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