Instruction Leaflet
Bedienungsanleitung
Hojas de instrucciones
Feuille d'instructions
Foglio d'istruzioni
Betjeningsvejledning
Instructies
Instruktionsfolder

Ultrasonic Cleaning Bath **GB**

Ultraschallreinigungsbad **D**

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
Kar til ultralydsrensning **DK**

Ultrasoonreiniger **NL**

Bad för ultraljudrengöring **SE**

1. Introduction

The key operation associated with this equipment is ultrasonic cleaning, so to give an understanding of the machine's operation, the principles of this technique are now described.

Ultrasonic cleaning depends upon the phenomena associated with acoustic cavitation. This is generated in a liquid medium as a direct result of subjecting it to continuous alternating pressure from piezo electric transducers. These generate acoustic vibrations in the frequency range 25 to 80 kHz.

Acoustic cavitation exists in two forms, stable and transient.

Stable cavities are bubbles of more or less uniform size, that are relatively permanent and will continue to oscillate for many cycles of applied pressure.

On the other hand transient cavities will generally exist for less than one cycle of pressure. During this cycle they will expand to many times their original size before collapsing violently.

Both cavities will contain either air or, more probably, vapour of the cleaning medium. Cleaning is effected by the gentle scrubbing action of the stable cavitation bubbles as well as the violent collapse of the transient cavitation bubbles.

Transient cavity pressure is dependent upon the frequency of the applied alternating pressure and hence, by selecting the ultrasonic transducer frequency, contamination ranging from the very heavy, as is found for example on engine components, to the very light, as found on semi-conductors or jewellery, can be removed.

In this way intricate objects having complex shapes can be cleaned without disassembly. The cavitation cleaning action occurs wherever the immersed parts come into contact with the cleaning medium in an operating ultrasonic tank.

Ultrasonic cleaning systems achieve extremely high standards of cleanliness in a fraction of the time required to obtain the same results by conventional methods such as hand cleaning, spray washing or similar mechanical means. They efficiently remove contaminants including polishing and lapping compounds, grinding abrasives, cutting oils and associated swarf, fluxes and general environmental contaminants from a wide range of components and sub-assemblies. Although the cavitation produced by ultrasonic waves is thorough, it will not damage delicate or intricate mechanisms.

Operation of an ultrasonic cleaning process requires a minimum of supervisory labour, and this may be unskilled. Multi-stage ultrasonic cleaning systems, designed for use with either organic solvents or aqueous solutions, may be readily incorporated into mass production lines. The development of automatic handling facilities for specific cleaning cycles is now well advanced.

2 Safety

2.1 Hazards

In normal operation, the main hazard to users of this equipment arises from possible exposure to excessive noise levels. There are subsidiary hazards associated with the handling of some of the more alkaline cleaning products and also from skin contact with hot cleaning liquids. Further significant hazards exist in respect of electrical failures and maintenance, where there is a risk of contact with live electrical circuits. The equipment is protected by a number of fuses, and in failure situations it is important that the cause is established and the malfunction rectified by a competent, qualified electrician.

2.2 Precautions

Observance of the following points is essential for the safe operation of the equipment.

- Ensure that the operating and maintenance instructions (see sections 3 and 4) are thoroughly understood before using the unit.
- The unit must be correctly wired and earthed, preferably using an earth leakage circuit breaker (ELCB). It must be set up for the correct supply voltage at its place of use.
- The cleaning bath must not be operated with liquid levels below the minimum level line (see Fig 1), and on no account must it be operated without liquid in the system.
- The unit must be disconnected from the mains before filling and emptying.

WARNING:

The warranty is invalidated if tank is operated with liquid level below indicated minimum level.

2.3 Cleaning Solutions

2.3.1 Type

It is to be emphasised that these units are only designed for use with neutral or mildly alkaline aqueous cleaning solutions. For Cleaning solutions, see the Cleaning/Contact treatment section of the current RS catalogue. The RS Technical Helpline will readily advise on the cleaning solution to be used for any particular application. The normal types of cleaners do not incorporate any recognised dangerous substances, but some of the more alkaline liquids require careful handling. Users must avail themselves of the manufacturers' technical and safety data sheets and ensure that the necessary precautions are taken in using the products.

For certain cleaning applications, the use of strongly alkaline solutions (containing caustic soda or sodium hydroxide) may be considered. Such cleaning media are compatible with the equipment, but where they are used the hazards to operatives are significantly increased, since caustic soda can cause severe skin burns and eye contact with this chemical can have very serious consequences.

As such, if strongly alkaline products are used, it is extremely important that the precautions detailed in the manufacturer's data sheets are followed.

In circumstances where rust or other corrosion deposits have to be removed from metal surfaces, the use of acidic cleaning solutions may be considered. If such cleaning is contemplated, it is essential to refer the matter to the RS Technical Helpline, since there is a real need to confirm the compatibility of the product with the equipment, so that corrosion is avoided.

RS units must not be used with organic solvents of the fluorocarbon or chlorinated hydrocarbon type, eg trichloroethylene. These solvents volatilise in air to give atmospheres that are harmful to operators.

2.3.2 Use of Cleaning Solutions

When filling, topping up or using an ultrasonic bath, contact of the cleaning solution with the eyes and the skin must be avoided. Impervious gloves (eg of PVC, polythene or rubber) should be worn routinely. Goggles (to BS 2092) should be worn when handling liquid concentrates and at all times when there is a risk of splashing with the working solution. Suitable protective clothing must be worn; if contaminated, clothing should be removed immediately and cleaned before re-use. The cleaning solutions tend to have a degreasing effect on the skin leaving it in a dry condition.

2.4 Safe Practices

Experience has shown that observance of the following simple rules will assist in securing the safer operation of cleaning units.

- DO NOT place unprotected hands in the cleaning solution, since the skin degreasing effect is enhanced by ultrasonic cavitation. Use a suitable component carrier.
- DO NOT drop heavy objects into a cleaning bath. There is a serious risk from splashing with liquid entering the eyes. Where such a risk is perceived, goggles must be worn.
- DO wear protective clothing, goggles and gloves and ensure that eye wash bottles are readily available.

Users must establish safe operating procedures for the running of ultrasonic baths, so that in their individual routine schedules, operations are conducted safely, without undue risk to personnel.

This can be achieved from a thorough understanding of the equipment and cleaning solvents and a disciplined application of the knowledge by managers, supervisors and operators to the particular cleaning situations. It is to be noted that despite the increasingly stringent standards imposed by safety legislation, there are many hundreds of RS Ultrasonic cleaning units in safe everyday use.

2.5 Noise Emission

Use of the ultrasonic cleaning system is attended by the emission of characteristic, high pitched noise consisting largely of harmonics and subharmonics of the basic working frequency (38 KHz). Many individual tests have shown that the total noise emitted from RS Ultrasonic baths is below the first action level of 85 (dBA) (Leq) (reference The Noise at Work Regulations 1989).

RS cleaning systems are configured such that the noise level from the equipment is below the statutory level under typical operating conditions. Despite this, there is a need for users to ensure that the noise intensity does not increase with time during routine use. In some circumstances, the provision of ear defenders for operators may be necessary. To combat excessive noise, much can be achieved by the sensible use of a tank lid when the equipment is on stand-by, since a good deal of noise is emitted through the liquid surface.

Particular care should be taken to avoid siting ultrasonic cleaning equipment near other noisy machinery, since the combined noise intensity could exceed the required safe limit. Where such a combination is unavoidable, then total noise level measurements must be made.

3. Operating Instructions

Strict compliance with the following operating procedures will ensure that the machine is utilised within its design parameters and will operate efficiently.

3.1 Siting and Installation

In order to ensure efficient operation of the system and to maintain a safe working environment for associated personnel, it is imperative that the following conditions are met:

- The unit must be placed on a level and secure base.
- Ensure that the unit is switched 'OFF' and that the heaters, if fitted, are 'OFF'.
- Check that the voltage supply is within the range on the serial number plate on the back panel
- Connect the input lead to a 13 amp plug and hence to the mains supply. Fuse rating to be as specified in Section 4.7. Connection must only be carried out by qualified electrical personnel.

The input cables must be connected to the plug as follows:

- Brown lead to the live terminal
- Blue lead to the neutral terminal
- Green/yellow lead to the earth terminal

It is very important that a good earth connection is made to a recognised earthing point.

When all these conditions are satisfied, the unit is ready for filling.

3.2 Filling of Unit

The appropriate safety precautions must be taken when filling the bath with cleaning solution (See section 2 - Safety).

As previously stated, RS ultrasonic tanks are designed for use with neutral or mildly aqueous solutions. These must be diluted with water to the concentrations recommended by the manufacturer. Occasional topping up is necessary to replace losses due to 'drag-out' of cleaning solution.

The level of cleaning liquid in the tank should be such that when the objects being cleaned are placed in the liquid, their displacement is not sufficient to cause an overflow.

On no account must the liquid depth be allowed to fall below the minimum level line inside the tank (see Fig. 1), or serious damage may be caused to the generator, transducers and heaters (where fitted).

To avoid thermal shock to the transducers, never pour hot liquids into a cold tank or vice versa. When filling a hot or cold tank, fill it approximately one quarter full with liquid of a similar temperature to that of the tank and complete the filling using either hot or cold liquid.

3.3 Start Up

- a) Ensure that all the conditions detailed in sections 3.1 and 3.2 are met.
- b) Switch the power on and turn the timer knob to set the processing time between 0 and 30 minutes. When processing is complete, the timer will switch the ultrasonic generators off.

3.4 Operation

Where the items to be cleaned are extremely soiled, they should be pre-washed or sprayed in a separate tank to remove bulk soil. An aqueous cleaning solution should be used and this may be the same as that used in the ultrasonic cleaning tank.

If the items to be cleaned are small and a number have to be done at one time, it is strongly advised that they should be immersed using baskets. The baskets must have perforated sides to allow the free movement of the cleaning liquid. An excessively deep load of items can cause masking in the ultrasonic tank and result in a slow or poor cleaning action. A shallow load will clean more rapidly and effectively.

The basket should be placed in the tank for a cleaning period that will normally be of 30 seconds to two minutes' duration, depending on the configuration of the items being cleaned and the degree of soiling. The basket must be withdrawn slowly to prevent excessive 'drag-out' of liquid.

The baskets supplied by RS components are specifically designed to prevent 'drag-out' of liquid and any similar work carrier devised by the user should achieve the same effect. Failure to allow adequate time for drainage, or removal of the basket at too great a speed, will result in excessive consumption of cleaning liquid.

On no account must the liquid level be allowed to fall below the minimum level line (see Fig 1) or serious damage may be caused to the generator and transducers.

3.5 Shutdown

On completion of a work period (eg at the end of the day or shift), or when no further cleaning is necessary, the unit must be shut down. Switch off the ON/OFF switch and remove the plug from the mains socket. Where a lid is supplied it should be left in position.

4. Maintenance

Maintenance is essential to keep the equipment operating efficiently. Basically it falls into two categories:

- routine, preventative operations that are performed by operators, and
- one-off operations to correct plant failures; these are usually beyond the scope of the normal operators and require the services of a qualified electrician or a technician

In all cases where the cause of failure is uncertain, reference should be made to RS Technical helpline.

The following schedules list the tasks that must be carried out, but the frequency of their implementation will be determined by the working environment of the machine, eg the type of components being cleaned, the nature of the cleaning liquid, the amount and type of contamination etc. A suggested initial maintenance frequency is annotated by each item.

Experience gained during operation of the machine will quickly establish a programme of maintenance that suits the particular working practice.

4.1 Safety Precautions

When any maintenance work is carried out on the machine, strict compliance with the following safety precautions must be enforced.

- Only qualified electrical personnel are permitted to remove any panels from the machine that will expose live terminals.
- Before any electrical connections are disturbed, the supply must be removed from the machine by removing the mains plug from the supply socket.
- Whenever the ultrasonic bath contains cleaning liquid, care must be taken to prevent contact with the skin and eyes. The wearing of impervious gloves and eye protection is necessary.

4.2 Routine Preventative Maintenance

In addition to the tasks listed, the operator should be constantly aware of the running condition of the machine. Any deterioration in performance should be investigated immediately and, where possible, the cause identified and the fault rectified.

Daily

- Check liquid level is above the minimum level (see Fig 1) and, if required, add further cleaning solution.
- Check condition of the cleaning fluid. As a general rule, evidence of unacceptable solution contamination is manifested by the visible presence of solid debris in the liquid and/or discolouration. Contaminated solution must be drained and replaced.

Monthly

- Check the security of the earth connection.
- Ensure that the mains plug is correctly mated in its socket.
- Visually examine electrical connections and supply leads particularly for looseness, fraying or signs of damage.

4.3 Draining and Cleaning the Plant


As work is being processed, cleaning results may begin to deteriorate because of the progressive build up of solid contamination: see Section 4.2, Daily maintenance (item b). Once this stage is reached the plant will need to be drained and cleaned. The frequency of cleaning will depend upon the work throughput and with experience can become a routine in particular applications.

4.5 Trouble Shooting Guide

This section is intended to be an aid in dealing with problems that can arise when using this equipment.

Effect	Cause	Action
High cleaning liquid consumption	Excessive liquid loss (drag-out)	- Use alternative basket design - Realign components in basket - Withdraw basket more slowly
Localised pitting corrosion	Heavy build-up of solid contamination including deposits	Drain and clean tank
Unsatisfactory cleaning	Solution very contaminated Masking by too many components	Change cleaning liquid Reduce number of components in the basket
Unit not working	Power failure	Check mains supply and all associated fuses

4.6 Replacement of Fuses

- Disconnect unit from power supply socket.
- Remove IEC power plug from rear of unit.
- Place medium flat-blade screwdriver into fuse holder slot then push in barrel and turn, then release and holder will spring outwards.
- Replace with correct fuses: 20mm ceramic quick acting fuses - (rating 220-240V) as follows:-

KC2	-	1A
KC3	-	1A
KC6	-	3.15A
KC14	-	5A
KC22	-	10A
- Push barrel back in and turn to lock in place, replace IEC plug.

As with filling the unit, great care should be taken when draining the plant to avoid eye and skin contact with liquid through splashing, spillage etc. During draining observe all safety features described in Section 2.

To drain the equipment proceed as follows:

- Switch off the plant and allow cleaning solution to cool.
- Drain the tank. The contaminated cleaning solution should be disposed of to a drain, diluting with copious volumes of water: (see Section 4.4 below).
- Thoroughly rinse the tank with water and empty the washings.
- Clean out the tank, paying particular attention to the base and corners. Any deposit or surface dullness will appreciably decrease the ultrasonic activity. If a cleaning agent is used, ensure that all traces are removed.
- The unit may now be refilled and brought back into operation. Ensure that the liquid level is above the minimum level (see Fig. 1).

4.4 Contaminated Liquid Disposal

The normal range of neutral and alkaline cleaners used in RS tanks are readily biodegradable and may be disposed of with simultaneous dilution to the public sewer. However, strongly alkaline (caustic soda) or strongly acidic solutions should be disposed of through an approved waste disposal contractor to a licensed site. These liquids should be stored in labelled plastic containers, although metal containers are suitable for caustic soda solutions.

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