USB/WiFi Streamer

Instruction Manual

Contents

- Features.....1
- Package Contents.....2
- Product overview and function.....3
- Led indicators.....4
- How to connect the USB/WiFi Streamer.....5
- Download and installation App.....6
- How to use “MicFi” app.....10

• INTRODUCTION

Wi-Fi Digital box streamer is a very innovative device that allows to stream from USB to Wi-Fi all kind of images taken from a digital USB device. It can be used with digital USB microscopes (i.e. Dino-Lite) but even with other devices like USB endoscopes or USB eyes. It breaks the conception of traditional USB digital microscope to realize the following functions: measurement, conservation, copy and transfer of images and video which are difficult for a traditional microscope.

You can use Wi-Fi Digital box even as power bank to charge your smartphone

• Features

Operation system	Windows, OS X, iOS, Android
Protocol	IEEE 802.11 b/g/n standard compliants
Distanza di Trasmissione Wireless	Up to 5 m
Transmit power	11n HT40 MCS7: +13 dBm, 11b CCK: +18 dBm, 11g OFDM: +15 dBm
Data rate	802.11n: up to 150Mbps , 802.11b: 1, 2, 5,5, 11Mbps , 802.11g: 6, 9, 12, 18, 24, 36, 48, 54Mbps
Frequency range	2.400 ~ 2.4835GHz
Power consumption	Max 2.5W
Li-ion battery	Continuous Working Time—Appx.3 Hours, Full—Appx.3 Hours
Power source	DC5.0V/1A

•DINO-LITE COMPATIBILITY

Series	Resolution	Led
Premier	1.3 Mpixels	NO
Premier R4	1.3 Mpixels	YES
Premier	5 Mpixels	YES
Edge	1.3 Mpixels SXGA	YES

Dino-Lite is a registered trade mark from third parties and is mentioned for descriptive use only

• PACKAGE CONTENTS

USB/WiFi Streamer


Power source (5V, DC 1A)

USB cable

• PRODUCT OVERVIEW & FUNCTION

- 1. Mini USB port** – This allows to be connected a mini USB cable carrying power from power source.
- 2. Power Switch** – This power switch turns On and OFF the USB/WiFi Streamer.
- 3. WiFi Indicator** – Read the information in the “LED indicators” section of this manual.
- 4. Charge Indicator** – Read the information in the “LED indicators” section of this manual.
- 5. Battery Indicator** – Read the information in the “LED indicators” section of this manual.
- 6. USB Interface** – This allows to connect USB microscopes, endoscopes and others USB devices

3

• LED INDICATOR

	LED COLOR	STATUS	DESCRIPTION
Charge Indicator	NO light	OFF	The device is powered off or otherwise is not getting power
Charge Indicator	GREEN	SOLID	The battery is charging
WiFi Indicator	NO light	OFF	The device is powered off
WiFi Indicator	BLUE	BLINKING	The device is transmitting or receiving data
WiFi Indicator	BLUE	SOLID	The device is powered on but is not transmitting or receiving data
Battery Indicator	NO light	OFF	The device is powered off or otherwise is not getting power
Battery Indicator	ORANGE	BLINKING	The battery needs urgent charge
Battery Indicator	ORANGE	SOLID	The device is powered on or it is charging

• HOW TO CONNECT USB/WIFI STREAMER

1. Connect the mini USB cable to USB/WiFi Streamer and to power source.

PLEASE MAKE A CHARGE OF AT LEAST 3HOURS.

DO NOT USE USB/WIFI STRAMER during the charge of the battery

2. Disconnect USB cable for the charge the battery and connect your USB device to USB interface (position n.6) and **power ON the USB/WiFi Streamer by Power switch (position n.2)**

• Download the app and install it

1. Install:
 - a) For **iPad** and **iPhone**: Search “**Mic-Fi**” in **App Store** to download and install it.
 - b) For **Smart-Phone** and **Tablet Android**: Search “**Mic-Fi**” in **Google Play** to download and install it.
 - c) For **Windows** and **Mac PC**: Download the Mic-Fi Software by www.mic-fi.com
2. After powered ON the USB/WIFI Streamer, wait for ~20 seconds until the led of Wi-Fi indicator  are solid.
3. Check the wireless connection on your Phone/Tablet or PC and a list of available networks will appear. Select the SSID “Wi-Fi2SCOPE_xxxxxx” network from the list to connect; the network “Wi-Fi2SCOPE_xxxxxx” requires a password: enter the password “**12345678**” and click Ok.
4. Once connected, please wait for about 5 seconds to open the “Mic-Fi” application to see your images. If there is no video display, please close and restart “Mic-Fi” app again. The Wi-Fi indicator LED  flashes and video stream signal should be displayed in the app. Video streaming can be broadcasted to many devices, but the frame rate of live video will decrease as the number of receiving devices increase.


Setting for Internet and Microscope’s Wi-Fi Connection at the same time

You can use the Wi-Fi network of the microscope and that of the PC / SmartPhone / Tablet at the same time

1. To access the internet by Wi-Fi connection at the same time the Wi-Fi microscope is working, open an internet browser, input <http://10.10.1.1> and click “enter”. A pop-up window will require you to enter a username and password, please enter: **admin** and **admin**, and then this window should appear:

2. Click  in the window to pop-up the screen:

Select "Wireless AP Client Mode", click  to pop-up the screen:

Click  till your destre network SSID appears, and click it to pop-up the screen:

Insert the password and click "Next". Now you can use USB/WIFI streamer and Internet at the same time

7

WiFi: Security setup

Open an internet browser, input <http://10.10.1.1> and click "enter". A pop-up window will require you to enter a username and password, please enter: **admin** and **admin**, and then this window should appear:

Click on

Click “Network Security Settings” to pop-up the screen on the right. Enter your desired Network Name (SSID) and password in “WPA-PSK Key” and “Re-Enter”, and click “Next”. The setting is done, now you can re-connect the new SSID Wi-Fi by entering the new password.

“MicFi” UTILIZZO DELL’APPLICAZIONE


Note: First to open the Mic-Fi app, connect the Wi-Fi microscope

“MicFi” on iPad/iPhone

Search “MicFi” in **App Store** to download and install it

Click the  icon to display video

Warm-up Tips: If the live video is black or slow, please restart “Mic-Fi” app again.

Capture photo or record video (The pixel resolution has to be set at 640x480 when recording video).

Videos and photos folder

Navigation icon - moving a dot left/up/down/right for calibration or measurement operation.

Measurement and calibration icon

Calibration;

Line;

Angle;



Circle (three dots)

Switch between photo taking and video recording

Resolution Setting. (When changing the resolution, please wait until the illumination LED is off and lights up again – it takes about 5 seconds).

Click for power on Snapshot

Add new folder

Measurement & Calibration


To obtain the accurate measurement data for each magnification, please calibrate before any measurement function.

Use the enclosed calibration ruler as to-be-observed object. Bring the image into focus by turning the "Focus & Magnification adjustment wheel";

Click  to pop-up the screen on the left. Click "Calibration" to back the live video, click two dots (moving the dot by clicking    ) which draw a your desired line, then click on the live video to show the following screen:

Enter the value of length and click “ok”:


Now the calibration is completed. Click on screen again, you can make measurement features accurately on the live video at this same focus. (Don't adjust again the “Focus & Magnification wheel”). If the new photo object is not in focus, please only adjust the distance between the microscope and the object until the photo is clear. The measurement data can be saved inside the taken photo.

“MicFi” on Android Tablet/Phone

Note: First to open the Mic-Fi app, connect the Wi-Fi microscope
Search “MicFi” in **Google play** to download and install it

Click the  icon to display video

Warm-up Tips: If the live video is black or slow, please restart “Mic-Fi” app again.

Capture photo or record video (The pixel resolution has to be set at 640x480 when recording video).

Videos and photos folder

Navigation icon - moving a dot left/up/down/right for calibration or measurement operation.

Measurement and calibration icon

Calibration;

Line;


Angle;

Circle (three dots)

Resolution Setting. (When changing the resolution, please wait until the illumination LED is off and lights up again – it takes about 5 seconds).

Snapshot

Power on snapshot.

Rename(V)

Rename the photo first to save it.

Default

Default setting

Measurement & Calibration

See iPhone/iPad measurement and calibration information

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