



## HypaCover Low Adherent Dressings - 5cm x 5cm, Pack of 10

- Absorbs fluid without sticking to skin
- Individually wrapped
- Suitable for abrasions, cuts and wounds

Low adherent, absorbent wound dressing helps draw exudates away from the wound. The film surface of this dressing helps prevent the pad from sticking to the wound surface and minimises trauma for casualty upon removal.

Use dressing retention sheets or microporous tape to hold the dressing in place. HypaCover Low Adherent Absorbent Dressings are suitable for abrasions, superficial cuts and lightly exuding wounds- also suitable for burns treatment.

### Product Specifications

SKU	D4957PK10
GTIN	5037331019482
dimensions	5x5cm
dressing_size	5x5cm
brand	hypacover



